

Fi-G Fi-P

Fi-G enables the shortest heating and melting time within 15s

Fi-G: Time for 0°C-230°C (S)

Woodpecker	Korean B brand	Chinese C brand	Chinese D brand
15	24	51	56

Fi-P



0.2 second heating time

High-performance conversion system ensures the efficient cutting




0.35mm bendable tip

Thinner tips can penetrate deep into the narrow root canal, providing you a better treatment effect.



High precision temperature control

Constant temperature prevents scalding



Fi-G



Quick heating within **15 seconds**

reduces time for waiting and increases your working efficiency.



0.5mm injection needle

Easy to reach the apical 3mm of canal brings better filling.



Prevent gutta-percha leakage

360 rotatable needle

Optimized structure, 360° uniform rotation, effectively avoiding gutta-percha leakage.

